

Cross Reference Guide New Lab Balances (Europe)



turning science **into solutions**

Find the Balances That
Best Fits Your Customers Needs!

1

What Kind of
Lab Balance is Your
Customer Looking for?




2

Find the
New Models!

3

Cross Check Superior
Features of the New
Lab Balance Models.



CPA Series	Replacement Higher Performance Level	Replacement Lower Performance Level
CPA26P*	-	-
CPA26P-OCE*	-	-
CPA225D*	-	-
CPA225D-OCE*	-	-
CPA324S	MSE324S-000-DU	-
CPA324S-OCE	MSE324S-1CE-DU	-
CPA224S	MSE224S-100-DU	SECURA224-1S
CPA224S-OCE	MSE224S-1CE-DU	SECURA224-1CEU
CPA124S	MSE124S-100-DU	SECURA124-1S
CPA124S-OCE	MSE124S-1CE-DU	SECURA124-1CE
CPA64	MSE124S-100-DU	SECURA124-1S
CPA64-OCE	MSE124S-1CE-DU	SECURA124-1CEU
CPA1003S	MSE1203S-000-DE	-
CPA1003S-OCE	MSE1203S-0CE-DE	-
CPA623S	MSE623S-100-DE	SECURA613-1S
CPA623S-OCE	MSE623S-1CE-DE	SECURA613-1CEU
CPA423S	MSE623S-100-DE	SECURA513-1S
CPA423S-OCE	MSE623S-1CE-DE	SECURA513-1CEU
CPA323S	MSE323S-100-DE	SECURA313-1S
CPA323S-OCE	MSE323S-1CE-DE	SECURA313-1CEU
CPA223S	MSE323S-100-DE	SECURA213-1S
CPA223S-OCE	MSE323S-1CE-DE	SECURA213-1CEU
CPA6202S	MSE6202S-100-DO	SECURA6102-1S
CPA6202S-OCE	MSE6202S-1CE-DO	SECURA6102-1CEU
CPA4202S	MSE4202S-100-DO	SECURA5102-1S
CPA4202S-OCE	MSE4202S-1CE-DO	SECURA5102-1CEU
CPA3202S	MSE4202S-100-DO	SECURA3102-1S
CPA3202S-OCE	MSE4202S-1CE-DO	SECURA3102-1CEU
CPA2202S	MSE2202S-100-DO	SECURA2102-1S
CPA2202S-OCE	MSE2202S-1CE-DO	SECURA2102-1CEU
CPA34001S	MSE36201S-000-DO	-
CPA34001S-OCE	MSE36201S-0CE-DO	-
CPA16001S	MSE20201S-000-DO	-
CPA16001S-OCE	MSE20201S-0CE-DO	-
CPA12001S	MSE20201S-000-DO	-
CPA12001S-OCE	MSE20201S-0CE-DO	-
CPA10001	MSE12201S-000-DO	-
CPA10001-OCE	MSE12201S-0CE-DO	-
CPA8201	MSE8201S-000-DO	SECURA6102-1S
CPA8201-OCE	MSE8201S-0CE-DO	SECURA6102-1CEU
CPA5201	MSE5201S-100-DO	SECURA5102-1S
CPA5201-OCE	MSE5201S-1CE-DO	SECURA5102-1CEU
CPA2201-OCE	MSE5201S-1CE-DO	SECURA5102-1CEU
CPA34000	MSE36200S-000-DO	-
CPA34000-OCE	MSE36200S-0CE-DO	-

* Still available until further notice



ED Series

	Replacement Higher Performance Level	Replacement Lower Performance Level
ED224S	SECURA224-1S	QUINTIX224-1S
ED224S-OCE	SECURA224-1CEU	QUINTIX224-1CEU
ED124S	SECURA124-1S	QUINTIX124-1S
ED124S-OCE	SECURA124-1CEU	QUINTIX124-1CEU
ED623S	SECURA613-1S	QUINTIX613-1S
ED632S-CW	SECURA613-1S	QUINTIX613-1S
ED623S-OCE	SECURA613-1CEU	QUINTIX613-1CEU
ED423S	SECURA513-1S	QUINTIX513-1S
ED423S-DS	SECURA513-1S	QUINTIX513-1S
ED423S-CW	SECURA513-1S	QUINTIX513-1S
ED423S-OCE	SECURA513-1CEU	QUINTIX513-1CEU
ED323S	SECURA313-1S	QUINTIX313-1S
ED323S-DS	SECURA313-1S	QUINTIX313-1S
ED323S-CW	SECURA313-1S	QUINTIX313-1S
ED323S-OCE	SECURA313-1CEU	QUINTIX313-1CEU
ED153	SECURA213-1S	QUINTIX213-1S
ED153-DS	SECURA213-1S	QUINTIX213-1S
ED153-CW	SECURA213-1S	QUINTIX213-1S
ED153-OCE	SECURA213-1CEU	QUINTIX213-1CEU

ED6202S	SECURA6102-1S	QUINTIX6102-1S
ED6202S-CW	SECURA6102-1S	QUINTIX6102-1S
ED6202S-OCE	SECURA6102-1CEU	QUINTIX6102-1CEU
ED4202S	SECURA5102-1S	QUINTIX5102-1S
ED4202S-CW	SECURA5102-1S	QUINTIX5102-1S
ED4202S-OCE	SECURA5102-1CEU	QUINTIX5102-1CEU
ED3202S	SECURA3102-1S	QUINTIX3102-1S
ED3202S-CW	SECURA3102-1S	QUINTIX3102-1S
ED3202S-OCE	SECURA3102-1CEU	QUINTIX3102-1CEU
ED2202S	SECURA2102-1S	QUINTIX2102-1S
ED2202S-CW	SECURA2102-1S	QUINTIX2102-1S
ED2202S-OCE	SECURA2102-1CEU	QUINTIX2102-1CEU
ED822	SECURA1102-1S	QUINTIX1102-1S
ED822-CW	SECURA1102-1S	QUINTIX1102-1S
ED822-OCE	SECURA1102-1CEU	QUINTIX1102-1CEU
ED8201	SECURA6101-1S	QUINTIX6101-1S
ED8201-CW	SECURA6101-1S	QUINTIX6101-1S
ED8201-OCE	SECURA6101-1CEU	QUINTIX6101-1CEU
ED5201	SECURA5101-1S	QUINTIX5101-1S
ED5201-CW	SECURA5101-1S	QUINTIX5101-1S
ED5201-OCE	SECURA5101-1CEU	QUINTIX5101-1CEU
ED2201	SECURA2101-1S	QUINTIX2101-1S
ED2201-CW	SECURA2101-1S	QUINTIX2101-1S
ED2201-OCE	SECURA2101-1CEU	QUINTIX2101-1CEU



TE Series	Replacement Higher Performance Level	Replacement Lower Performance Level
TE214S	QUINTIX224-1S	PRACTUM224-1S
TE214S-OCE	QUINTIX224-1CEU	-
TE124S	QUINTIX124-1S	PRACTUM124-1S
TE124S-OCE	QUINTIX124-1CEU	-
TE64	QUINTIX64-1S	PRACTUM64-1S
TE64-OCE	QUINTIX64-1CEU	-
TE313S	QUINTIX313-1S	PRACTUM313-1S
TE313S-DS	QUINTIX313-1S	PRACTUM313-1S
TE153S	QUINTIX213-1S	PRACTUM213-1S
TE153S-DS	QUINTIX213-1S	PRACTUM213-1S
TE3102S	QUINTIX3102-1S	PRACTUM3102-1S
TE1502S	QUINTIX2102-1S	PRACTUM2102-1S
TE612	QUINTIX612-1S	PRACTUM612-1S
TE512-OCE	QUINTIX612-1CEU	-
TE412	QUINTIX412-1S	PRACTUM412-1S
TE212	QUINTIX412-1S	PRACTUM412-1S
TE6101	QUINTIX6101-1S	PRACTUM6101-1S
TE4101	QUINTIX5101-1S	PRACTUM5101-1S
TE2101	QUINTIX2101-1S	PRACTUM2101-1S
TE601	QUINTIX2101-1S	PRACTUM2101-1S
TE12000	QUINTIX6100-1S	PRACTUM6100-1S
TE6100	QUINTIX6100-1S	PRACTUM6100-1S
TE4100	QUINTIX5100-1S	PRACTUM5100-1S



M-Pact Series

	Replacement Higher Performance Level	Replacement Lower Performance Level
AX224	QUINTIX224-1S	PRACTUM224-1S
AX124	QUINTIX124-1S	PRACTUM124-1S
AX623	QUINTIX613-1S	PRACTUM513-1S
AX423	QUINTIX513-1S	PRACTUM513-1S
AX6202	QUINTIX6102-1S	-
AX4202	QUINTIX5102-1S	PRACTUM3102-1S
AX2202	QUINTIX2102-1S	PRACTUM2102-1S
AX822	QUINTIX1102-1S	PRACTUM1102-1S




M-Power Series

	Replacement Higher Performance Level	Replacement Lower Performance Level
AZ214	PRACTUM224-1S	ENTRIS224-1S
AZ124	PRACTUM124-1S	ENTRIS124-1S
AZ64	PRACTUM64-1S	ENTRIS64-1S
AZ313	PRACTUM313-1S	ENTRIS323-1S
AZ153	PRACTUM213-1S	ENTRIS153-1S
AZ3102	PRACTUM3102-1S	ENTRIS3202-1S
AZ1502	PRACTUM2102-1S	ENTRIS2202-1S
AZ612	PRACTUM612-1S	ENTRIS822-1S
AZ212	PRACTUM412-1S	ENTRIS822-1S
AZ6101	PRACTUM6101-1S	ENTRIS8201-1S
AZ4101	PRACTUM5101-1S	ENTRIS5201-1S
AZ2101	PRACTUM2101-1S	ENTRIS2101-1S



How to Find the Right Balance?

Sartorius' new lab balance models are designed customer specific and for their individual needs and applications. The below cross reference guide will help you choose the balance that directly fits to your needs.



Features	Cubis® MSE	Secura®	Quintix®	Practum®	Entris®
Modularity	•				
Password Protection	•	•			
Monitoring of the minimum sample weight acc. to USP, Chapter 41		•			
Sample & Batch Identifiers	•	•			
GLP GMP compliant Print-Out	•	•	•		
Direct Data Transfer into Windows® Applications		•	•	•	
Blocked Data Transfer		•			
Easy to Clean	•	•	•	•	•
Weighing Ranges					
Advanced (including µg)	•				
Standard	•	•	•	•	
Basic	•	•	•	•	•
Applications					
Advanced Lab Applications		•	•		
Standard Lab Application	•	•	•	•	
Basic Lab Applications	•	•	•	•	•
Leveling					
Automatic leveling	•				
Level monitoring and control	•	•			
Manual leveling	•	•	•	•	•
Interface					
Bluetooth®	•				
USB	•	•	•	•	
RS232	•	•*	•*	•*	•
Draft Shield					
Motorized Draft Shield	•				
Motorized Draft Shield with deionizer	•				
Manual Draft Shield	•	•	•	•	•
Operation					
Keys, Liquid crystal display (black white)	•				•
Colored touchscreen		•	•	•	
Calibration					
Internal Adjustment with isoCAL	•	•	•		
External Adjustment	•	•	•	•	•

* with using adapter cable

Sartorius Lab Instruments GmbH & Co. KG
Weender Landstrasse 94–108
37075 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius.com



◀▶ www.sartorius.com